

# PerOx (TOS/TOC) Kit

*Photometric test system for the determination of the total oxidative status/capacity (TOS/TOC) in EDTA plasma, serum and other biological samples*

Valid from 2019-01-01

**REF** KCR5100

$\Sigma$   
96

+2°C  
+8°C

-20°C  
CAL  
CTRL 1  
CTRL 2

**RUO**



**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)



# Table of Contents

<b>1. INTENDED USE</b>	<b>2</b>
<b>2. INTRODUCTION</b>	<b>2</b>
<b>3. MATERIAL SUPPLIED</b>	<b>2</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>2</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>3</b>
<i>Calibrator and controls</i>	3
<i>Storage of the other reagents</i>	3
<i>Preparation of the reaction buffer mixture</i>	3
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>4</b>
<i>EDTA plasma and serum</i>	4
<i>Cell culture</i>	4
<b>7. ASSAY PROCEDURE</b>	<b>4</b>
<i>Principle of the test</i>	4
<i>Test procedure</i>	5
<b>8. RESULTS</b>	<b>5</b>
<b>9. LIMITATIONS</b>	<b>6</b>
<b>10. QUALITY CONTROL</b>	<b>6</b>
<i>Reference ranges</i>	6
<i>Controls</i>	7
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>7</b>
<i>Precision and reproducibility</i>	7
<i>Detection limit</i>	7
<b>12. PRECAUTIONS</b>	<b>7</b>
<b>13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>8</b>
<b>14. REFERENCES</b>	<b>8</b>

## 1. INTENDED USE

This photometric Immundiagnostik assay is intended for the quantitative determination of the total oxidative status/capacity (TOS/TOC) in EDTA plasma, serum and cell culture supernatants. For research use only only. Not for use in diagnostic procedures.

## 2. INTRODUCTION

An overproduction of oxygen radicals or insufficient antioxidative capacity leads to a dangerous imbalance in the organism. This starts pathological mechanisms.

The **PerOx** assay is fast, reliable and easy to perform. **Total lipid peroxides** are measured. Because of a direct correlation between oxygen radicals and lipid peroxides, it is possible to measure and characterize the **oxidative status/oxidative stress** in biological fluids.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KCR5100	CAL	Calibrator (lyoph. 0.25 ml)	3 vials
KCR5100	CTRL1 CTRL2	Control 1 and 2 (lyoph. 0.25 ml)	3 vials each
KCR5100	REABUF A	Reaction buffer A	25 ml
KCR5100	REABUF B	Reaction buffer B	1 ml
KCR5100	ENZ	Enzyme solution	50 µl
KR0005.15	RECSOL	Reconstitution solution	15 ml
KCR5100	STOP	Stop solution	15 ml
KCR5100	PLATE	Microtiter plate	1 plate

Individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

For reorders of single components, use the catalogue number followed by the label as product number.

## 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Vortex
- Various pipettes
- Incubation chamber for 37°C
- Microtiter plate reader (required filters see chapter 7)

## 5. STORAGE AND PREPARATION OF REAGENTS

### *Calibrator and controls*

The **lyophilised calibrator** (CAL) and the **lyophilised controls 1 and 2** (CTRL1 and CTRL2) are stable at **-20°C** until the expiry date stated on the label. Before use, they have to be reconstituted with each **250 µl reconstitution solution** (RECSOL). Allow the vial content to dissolve for 5 min and then mix thoroughly by vortexing. Aliquots of the **calibrator** (reconstituted CAL) and the **controls 1 and 2** (reconstituted CTRL1 and CTRL2) can be stored at **-20°C** for **4 weeks**. **Avoid repeated thawing and freezing**. The concentration of calibrator and controls slightly changes from lot to lot. The exact concentration is stated on the label of CAL and in the specification of the controls, respectively.

### *Storage of the other reagents*

Reaction buffer B (REABUF B) has to be stored at **2–8°C in the dark**.

Test reagents are stable until the expiry date stated on the label when stored at **2–8°C**.

### *Preparation of the reaction buffer mixture*

To avoid losses, the enzyme solution (ENZ) should be centrifuged prior to use. After use, the vial has to be tightly closed to avoid contamination or evaporation (e.g. with parafilm).

The **reaction buffer mixture** must be prepared directly before use:

- 5 ml** reaction buffer A (REABUF A)
- + **100 µl** reaction buffer B (REABUF B)
- + **5 µl** enzyme solution (ENZ)

**Important:** The amounts mentioned above are sufficient for 40 tests. For other sample numbers, the buffer volumes must be adjusted accordingly.

**Important:** The reaction buffer mixture **cannot be stored**. The reaction buffers (REABUF A, REABUF B) and the enzyme solution (ENZ) are stable at 2–8°C until the expiry date stated on the label.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *EDTA plasma and serum*

Venous fasting blood is suited for this test system. EDTA plasma should be preferred because in serum, a time dependent increase in peroxide concentration is observed. EDTA plasma is stable at  $-20^{\circ}\text{C}$  for at least 4 weeks.

During serum preparation, it is important not to exceed 30 min at room temperature for clotting. Afterwards, the serum should be stored at  $-20^{\circ}\text{C}$  up to the measurement.

Samples with visible amounts of precipitates (mostly cryoproteins) should be centrifuged (at least 5 min at 10000g) prior to measurement. The supernatant is used in the test.

### *Cell culture*

In principle, it is possible to determine the PerOx concentration in cell culture supernatants. To test whether ingredients of the cell culture medium affect the measurement or not, we recommend the following approach:

- Dilute  $10\ \mu\text{l}$  of 30%  $\text{H}_2\text{O}_2$  concentrate with  $1000\ \mu\text{l}$  of reaction buffer A (REABUF A) =  $S_0$ .
- Dilute  $5\ \mu\text{l}$  of  $S_0$  with  $1000\ \mu\text{l}$  of reaction buffer A (REABUF A) =  $S_1$ .
- Preparation 1: Pipet  $10\ \mu\text{l}$  cell culture medium into one well.
- Preparation 2: Pipet  $10\ \mu\text{l}$  aqua bidest. into another well.
- Add each  $10\ \mu\text{l}$   $S_1$  to both preparations, then add each  $100\ \mu\text{l}$  reaction buffer A (REABUF A) and  $100\ \mu\text{l}$  reaction buffer mixture (see preparation above). Incubate for 5 min. Add  $50\ \mu\text{l}$  stop solution (STOP) and measure immediately at 450 nm in a microtiter plate reader.

### **Evaluation of the cell culture medium results**

A ratio of  $\text{OD}_{\text{preparation 1}} : \text{OD}_{\text{preparation 2}} > 0,8$  demonstrates that there are no major disturbing factors in the tested cell medium, and the assay can be performed.

## 7. ASSAY PROCEDURE

### *Principle of the test*

The determination of the peroxides is performed by the reaction of a peroxidase with peroxides in the sample followed by the conversion of TMB to a colored product.

After addition of a stop solution the samples are measured at 450 nm in a microtiter plate reader. The quantification is performed by the delivered calibrator.

### *Test procedure*

#### **Procedural notes**

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according the enclosed manual.

#### **Test procedure**

The microtiter plate is ready to use.

**Important:** To ensure the reproducibility of the measurement, the given incubation time and temperature must strictly be followed.

1.	Pipet <b>10 µl</b> of sample, calibrator (CAL), controls (CTRL1, CTRL2) and blank/reconstitution solution (RECSOL) in appropriate wells.
2.	Add <b>100 µl</b> of reaction buffer A (REABUF A).
3.	<b>Measurement 1:</b> Read the absorption of the samples in the ELISA reader at 450 nm.
4.	Add <b>100 µl</b> of reaction buffer mixture.
5.	Incubate for <b>15 min</b> at <b>37 °C</b> .
6.	Add <b>50 µl</b> stop solution (STOP).
7.	<b>Measurement 2</b> is performed immediately after addition of the stop solution (STOP) at 450 nm in the ELISA reader.

## **8. RESULTS**

The difference between measurement 1 and 2 is directly proportional to the peroxide content of the sample:

For evaluation, the estimated OD values of the first measurement are subtracted from the optical densities of the second measurement to obtain the  $\Delta OD$  values of sample, calibrator, controls and blank.


The concentrations of samples and controls are calculated using the calibrator (see label for concentration) and the following formula:

$$\text{Sample conc. } [\mu\text{mol/l}] = \frac{\Delta OD \text{ sample} - \Delta OD \text{ blank}}{\Delta OD \text{ calibrator} - \Delta OD \text{ blank}} \times \text{conc. calibrator } [\mu\text{mol/l}]^*$$

\* see label for concentration

A prepared Excel evaluation file can be requested by Immundiagnostik AG.

The following linear standard curve is for demonstrational purpose only. Because of the linearity in the chosen concentration range, one-point calibration using the included calibrator is sufficient.



## 9. LIMITATIONS

Whole blood cannot be used.

Strong haemolytic and lipaemic samples often show pathological concentrations. We do not recommend analysis of such samples.

The use of heparin plasma results in wrong results. Therefore, heparin plasma cannot be used in this assay.

## 10. QUALITY CONTROL

### *Reference ranges*

We recommend each laboratory to establish its own reference ranges. The values mentioned above are only for orientation and can deviate from other published data.



## Controls

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-assay CV**

2.94 % (162  $\mu\text{mol/l}$ ) [n = 6]

#### **Inter-assay CV**

6.63 % (136  $\mu\text{mol/l}$ ) [n = 10]

6.85 % (389  $\mu\text{mol/l}$ ) [n = 10]

### *Detection limit*

7  $\mu\text{mol/l}$

## 12. PRECAUTIONS

- Calibrator (CAL) and controls (CTRL) are based on human plasma. They were tested and found to be negative for HIV and hepatitis B. However, for safety reasons, all kit components should be treated as potentially infectious.
- The stop solution consists of diluted sulphuric acid ( $\text{H}_2\text{SO}_4$ ), a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- Reagents should not be used beyond the expiration date stated on the kit label.












### 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the kit package are for research use only.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for laboratories should be followed.
- Control samples should be analyzed with each run.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

### 14. REFERENCES

1. Hildebrandt, W. et al., 2002. Effect of N-acetyl-cysteine on the hypoxic ventilatory response and erythropoietin production: linkage between plasma thiol redox state and O(2) chemosensitivity. *Blood*, **99**(5), pp.1552–5.
2. Reichenbach, J. et al., 2002. Elevated oxidative stress in patients with ataxia telangiectasia. *Antioxidants & redox signaling*, **4**(3), pp.465–9.
3. Schimke, I. et al., 2001. Decreased oxidative stress in patients with idiopathic dilated cardiomyopathy one year after immunoglobulin adsorption. *Journal of the American College of Cardiology*, (1), pp.178–83.

**Used symbols:**

	Temperature limitation		Catalogue Number
	For research use only		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		